


Botany 2014 Next Generation Sequencing Workshop


1. Sequencing Technology
2. Genome Reduction
3. Sequence Processing
4. Bionformatics
5. What you really need to know.

Improvements in the rate of DNA sequencing over the past 30 years




MR Stratton et al. (2009) The cancer genome. Nature **458**, 719-724

Cost per Raw Megabase of DNA Sequence



Developments in High Throughput Sequencing




© Lex Nederbragt 2013

<http://flxlexblog.wordpress.com/2013/10/01/developments-in-next-generation-sequencing-october-2013-edition/>


Decoding long nanopore sequencing reads of natural DNA

Andrew H Laszlo¹, Ian M Derrington¹, Brian C Ross¹, Henry Brinkerhoff¹, Andrew Adey², Ian C Nova¹, Jonathan M Craig¹, Kyle W Langford¹, Jenny Mae Samson¹, Riza Daza², Kenji Doering¹, Jay Shendure² & Jens H Gundlach¹



Each of the 256 possible quadromers has a unique current value

In Sequence 2013 Survey: Illumina Pulls Further Ahead




<http://www.genomeweb.com/sequencing/sequence-2013-survey-illumina-pulls-further-ahead-interest-oxford-nanopore-remain>

Illumina Sequencing


released in 2007 by Solexa

Template Type	Sequencing Method	Imaging Method
Clonally amplified by solid phase amplification	Sequencing by synthesis with cyclic reversible termination	Four color imaging of single events using fluorescence



http://www.illumina.com/systems/hiseq_2000.ilmn

Illumina Sequencing



input:

- 1 ng – 1ug DNA
- 1 ug total RNA

“Classic” Illumina Library Prep

genomic DNA, RNA
BACs, amplicons ...

DNA fragmentation


end-repair

A-tailing

adapter-ligation


size selection

enrichment



Illumina Sequencing

1 lane = 200 million reads



(a) Flow cell channels reagents between glass plates with millions of anchored DNA fragments for sequencing.

(b) Cross section of flow cell showing second strand synthesis of anchored template strands. Laser light activates base-specific emission colour from modified nucleotides incorporated at the end of each cycle. (c) Imaging the surface of the flow cell after each cycle records the base added at each 'spot' of DNA. (d) Image analysis 'reads out' the sequence of incorporated second strand bases at each spot.

CLOSE X

Illumina Sequencing

HiSeq

Run time	6-12 days
Read length	125+125 bp
Yield/lane	60 Gbp

MiSeq

Run time	1-3 days
Read length	300+300 bp
Yield/lane	10-15 Gbp

New Illumina Sequencing Platforms Announced January, 2014

“Factory” Scale

HiSeq X

1.8 Tb/ 3 days

150 bp paired ends

Patterned Flow Cells (billions of nanowells)

Only Human Genome Sequencing Supported!

Sold in Sets of 10...



Desktop Model

NextSeq 500

up to 400 million reads

Two-Channel Base Detection (faster run time)

75 or 150 bp single or paired ends



NGS Error Rates

2014 NGS Field Guide. Travis Glenn. www.molecularrecologist.com

Platform	Primary Errors	Single-pass Error Rate (%)	Final Error Rate (%)	Notes
3730xl (capillary)	Substitution	0.1-1	0.1-1	
454	Indel	1	1	
Illumina	Substitution	~0.1	~0.1	≥ 75-85% of reads
SOLiD	A-T bias	~5	≤0.1	2x-3x sequencing
Ion Torrent	Indel	~1	~1	0.46-2.4%
PacBio RS	Indel	~13	≤1	consensus of 3 reads
Oxford Nanopore	Deletion	≥4	4	press release only

How to Deal with Overcapacity


Multiplexing

Addition of a unique sequence identifier (barcode or index) allowing multiple samples to be run together on a single flow cell lane.

Internal indexes (Cronn et al 2008)

External indexes (extra round of sequencing)


- Illumina 24 single; 96 dual
- NEB 24 single, 96 dual
- Nextflex 48 single, 192 dual
- NuGen 384 dual



Multiplexing Example


MiSeq v.3 November, 2013
20.3 million paired end reads
69 *Fragaria virginiana* BACs

Ideal Equimolar Results



average = 295,000

Actual Results



1200-13.8 million reads
median = 236,000

Outsourcing Library Prep

1. Fairly standardized (although diverse library options exist)
2. Expensive equipment (sonication, quantification)
3. Automation leading to lower prices



Botany 2014 Next Generation Sequencing Workshop

1. Sequencing Technology
2. Genome Reduction
3. Sequence Processing
4. Bionformatics
5. What you really need to know.

Whole Genome Sequencing (not covered here)

Genome Reduction

- A. PCR, fosmids, BACs
- B. Low Coverage (=Genome Skimming)
- C. Restriction Digest Methods (GBS, RAD, etc.)
- D. Target Capture (Hyb-Seq)
- E. Transcriptome (RNA-Seq)


Cronn et al. 2012. Targeted enrichment strategies for next-generation plant biology. *American Journal of Botany* 99:291-311.

Straub et al. 2012. Navigating the tip of the genomic iceberg: next-generation sequencing for plant systematics. *American Journal of Botany* 99:349-364.

Lemmon & Lemmon. 2013. High-throughput genomic data in systematics and phylogenetics. *Annual Review of Ecology, Evolution, and Systematics* 44:99-121.

Ellegren H. 2014. Genome sequencing and population genomics in non-model organisms. *Trends Ecol. Evol.* 29:51–63.

Genome Reduction Approaches



A. Amplicons and BACs

= genome reduction prior to NGS library preparation

1. Amplicons Fluidigm automation

Uribe-Convers et al 2014. A Long PCR–Based Approach...
Appl. Plant Sci. 2:1300063.

48 amplicons X 48 samples = \$600


1000+ can be sequenced on on MiSeq run = \$1500

2. Fosmids (10-40 kbp) and BACs (100-150 kbp)

Relatively expensive (\$10,000 for a BAC library)

Standard approach in large (well-funded) genome projects.


B. Genome Skimming




Straub et al. 2012. American Journal of Botany 99:349-364.

Genome Assembly

Reference guided assembly



De novo assembly



Genome Coverage




Low Coverage




High Coverage

Plant Cells




©Kendall-Hunt Publishing Company

Plant Genomes



Arabidopsis thaliana




Chloroplast
genome 154
kbp

Mitochondrial
genome
367 kbp

Arabidopsis Genome
Initiative. 2000. *Nature*
408:796-815.


Genome Skimming



Straub et al. 2012. American Journal of Botany 99:349-364.


Genome Skimming

Important:
Don't use more
sequence than needed.
100-125x is all that is
needed for organelles



C. Restriction Digest Approaches (RAD-Seq, GBS, 2b-RAD)

<http://www.maizegenetics.net/>



Buckler Lab for Maize Genetics and Diversity

A USDA-ARS Lab with Cornell's Institute for Genomic Diversity

Main Menu

- [Home](#)
- [People](#)
- [Research](#)
- [Bioinformatics](#)
- [Publications](#)
- [Germplasm](#)
- [Jobs](#)
- [About Us](#)
- [Links](#)
- [Contact Us](#)
- [Site Map](#)
- [Repository](#)

Search Site

GBS Overview

Genotyping by sequencing (GBS) is a simple highly-multiplexed system for constructing reduced representation libraries for the Illumina next-generation sequencing platform developed in the Buckler lab by [Rob Elshire](#). Key components of this system are: reduced sample handling, fewer PCR and purification steps, no size fractionation and inexpensive barcoding. We use restriction enzymes to reduce genome complexity and avoid the repetitive fraction of the genome.

Quick Links:

- [GBS Bioinformatics](#)
- [Workshop Videos](#)
- [FAQ](#)
- [GBS Method Paper](#)
- [Presentation on GBS](#)
- [96 Plex GBS Protocol](#)
- [Dilution Calculator](#)
- [Bar Coded Adapter Generator](#) (outside link)
- [384 Plex ApeKI Adapters](#) (Updated May 11, 2012 to correct two bad bar codes.)

Advantages of Restriction Digest approaches

1. Can obtain thousands of SNPs without a reference genome
2. Inexpensive library prep = hundreds to thousands of individuals
3. Discovery of candidate loci associated with phenotypic traits (quantitative trait loci = QTL) in natural populations (genome wide association studies = GWAS)
4. Analytical pipelines available
5. Hundreds of studies published to date.

Why I Don't Like Restriction Digest approaches

1. Generally requires at least 1 μg of good quality DNA.
2. Restriction digestion adds another variable to the library prep.
 - a. Restriction site polymorphism results in missing data.
 - b. Short fragments are cut less efficiently.
3. Does not target specific genes or SNPs.
4. References are easily obtained, and soon to be widely available.
5. Loci are not transferrable among species.
6. High potential for ascertainment bias.

Restriction Site Bias





Fig. 8 Results for the *in silico* digests 102 *Drosophila melanogaster* genomes. (A) Proportion of sites located in distinct regions of the genome when *in silico* digests are performed with different enzyme recognition sequences. GC-rich recognition sequences sample more exons, whereas AT-rich recognition sequences sample comparatively more introns and intergenic regions. 'Random' values are calculated from fragments selected at random throughout the genome. (B) Box plots of true π for regions sampled by enzymes with different recognition sequences. (C) The median true π (solid line) and estimated π (dashed line) as a function of chromosome sampling depth for three different recognition sequences. (D) Median of the ratio of estimated π to true π as a function of the number of sampled chromosomes. Dark blue, purple and cyan lines represent the three different restriction enzymes used in the *in silico* digest of the *D. melanogaster* genomes, and the dotted black line is from simulations with $\rho = 0.1$ per bp, $\theta = 0.01$ per bp

D. Targeted Sequence Capture = Hyb-Seq



12-24 individuals per hybridization

0.45-2 Mbp / individual

48 individuals per MiSeq lane

=250-1000 nuclear genes
or 6500-20,000 SNPs

384 individuals per HiSeq lane

~50% efficiency
= 50% off-target genome
skimming

\$50-\$75 per individual

Weitemier et al. 2014 APPS, in press.

Advantages of Hyb-Seq

1. A single laboratory procedure and bioinformatics pipeline can be used for phylogenetics (deep and shallow), population genetics and genetic linkage mapping.
2. A relatively distant (e.g. plant family) reference can be used.
3. Candidate genes can be targeted.
4. Can be scaled from 250 genes to an entire exome (25,000 - 30,000 genes).
5. Data sets can be easily combined and extended (in contrast to GBS, SNP chips).
6. Minor potential for ascertainment bias.

Why not sequence the same genes?

BMC Genomics



This Provisional PDF corresponds to the article as it appeared upon acceptance. Fully formatted PDF and full text (HTML) versions will be made available soon.

Identification, characterization, and utilization of single copy genes in 29 angiosperm genomes

BMC Genomics 2014, **15**:504 doi:10.1186/1471-2164-15-504

Fengming Han (hanfengming2004@163.com)

Yong Peng (ypeng@implad.ac.cn)

Lijia Xu (ljxu@implad.ac.cn)

Peigen Xiao (xiaopg@public.bta.net.cn)

3000-11000 (8-35%) single copy genes in 29 angiosperms

12 single copy genes shared across all 29 species

Genome Reduction

- A. PCR, fosmids, BACs
- B. Low Coverage (=Genome Skimming)
- C. Restriction Digest Methods (GBS, RAD, etc.)
- D. Target Capture (Hyb-Seq)
- E. Transcriptome (RNA-Seq)

RNA-SEQ

For Everyone

<http://rnaseq.uoregon.edu/index.html>

Design
Experiment

1. Carefully design
the experiment.

Purify
RNA

2. Isolate and purify
input RNA.

Prepare
Libraries

3. Convert the RNA
to cDNA and add
sequencing adapters.

Sequence

4. Sequence cDNAs
using one of the
available NGS platforms.

Analysis


5. Analyze the resulting
short-read sequences.

transcriptome reference
genome reduction





de novo assembly


quantify gene expression




read mapping



Experimental Design



RNA-Seq Technical vs. Biological Replicates




Busby et al. 2011. Expression divergence measured by transcriptome sequencing of four yeast species. BMC Genomics 12: 635.

Scotty - Power Analysis for RNA Seq Experiments


Scotty is a tool to assist in the designing of RNA Seq experiments that have adequate power to detect differential expression at the level required to achieve experimental aims.

<http://bioinformatics.bc.edu/marthlab/scotty/scotty.php>



These have too little power and too much measurement bias

Count Normalization




RPKM = reads per kilobase per million reads (most popular method, performs the worst)

All of these methods assume a constant denominator (total RNA per cell)

Dillies et al. 2013. A comprehensive evaluation of normalization methods for Illumina high-throughput RNA sequencing data analysis. *Briefings in Bioinformatics* 14: 671-683.

Total RNA Amount Varies Among Samples

different developmental stages
different ploidy levels
(Ilut et al, 2012)




Potential Solutions:


Scale by Total RNA / Cell

Include Internal RNA Standards

Aanes et al. 2014. Normalization of RNA-Sequencing Data from Samples with Varying mRNA Levels. PLOS ONE 9: e89158.



RNA and proteins not perfectly correlated



1 = leaf base
 4 = 3-4 cm from base
 9 = 9-10 cm from base
 14 = leaf tip (13-14 cm)

Ponnala et al 2014. The Plant Journal 78:424-440.

RNA-Seq

Pros:

Powerful method with multiple applications:

SNPs, differential expression, genome annotation, phylogenomics, alternative splicing, RNA editing

Cons:

Requires living tissue

RNA extraction is more involved than DNA

RNA expression biases (biological and analytical) are far from understood

Analysis of isoforms and paralogs is challenging



Botany 2014 Next Generation Sequencing Workshop


1. Sequencing Technology
2. Genome Reduction
3. Sequence Processing
4. Bionformatics
5. What you really need to know.

Sequence Processing

Potential Contaminants and Biases

1. Low Quality Bases
2. Adapters
3. Duplicate Reads
4. Uneven Sequencing Depth
5. Fragment Size Variation
6. Read Overlap
7. Biological Contamination





Sequence Processing

Quality Assessment Tools:

Illumina BaseSpace <https://basespace.illumina.com/runs>

FastQC

<http://www.bioinformatics.babraham.ac.uk/projects/fastqc/>

SGA Simpson, 2014. Bioinformatics
 uses kmers to estimate:
 paired end insert sizes
 heterozygosity
 repeat content
 genome size

Sequence Processing

A Popular Read Trimming & Adapter Removal Tool:

Trimmomatic (Bolger et al 2014)

<http://www.usadellab.org/cms/?page=trimmomatic>

- trims bases at start & end of reads, using a threshold quality score
- also can use a sliding window average of base quality to trim
- filters reads below a minimum length after trimming
- keeps track of orphaned pairs
- short adapters can be trimmed in palindrome mode

Phred = 20 (1% error rate) is commonly used

But is this too stringent? (MacManes 2014)

Sequence Processing

Duplicate Removal

Optical duplicates = split clusters

PCR duplicates from final “enrichment” step of library prep

Best to limit to absolute minimum of cycles (6-12)

Deep sequencing (>500X) will result in “natural duplication”

Should these be removed? Differences of opinion exist.

See Zhou et al. Bioinformatics January 2, 2014.

Yes for SNPs; No for counting [but incorporate into model]

Tools for Duplicate Read Removal:

Picard (picard.sourceforge.net; requires mapping)

FastUniq (Xu et al 2012; requires paired sequences)

fastq_collapse.py (Weitemier, unpublished)

Sequence Processing

k -mer = substrings of length k in DNA sequence data

ACGTCACG

acg tca


cgt cac five 3-mers

gtc **acg**

Counting the occurrences of all such substrings in a sequence generates a k -mer frequency distribution:

(reverse complements counted as a single k -mer)

Melsted & Pritchard (2011) Efficient counting of k -mers in DNA sequences using a bloom filter. BMC Bioinformatics 12:333.



Gray = errors and SNPs

Sequence Processing

Digital Normalization

Brown et al (2014)

1 error => up to k erroneous k -mers

AAAAAGAAAA

5-mer example

Table 1. Digital normalization to $C=20$ removes many erroneous k -mers from sequencing data sets. Numbers in parentheses indicate number of true k -mers lost at each step, based on reference.

Data set	True 20-mers	20-mers in reads	20-mers at $C=20$	% reads kept
Simulated genome	399,981	8,162,813	3,052,007 (-2)	19%
Simulated mRNAseq	48,100	2,466,638 (-88)	1,087,916 (-9)	4.1%
<i>E. coli</i> genome	4,542,150	175,627,381 (-152)	90,844,428 (-5)	11%
Yeast mRNAseq	10,631,882	224,847,659 (-683)	10,625,416 (-6,469)	9.3%
Mouse mRNAseq	43,830,642	709,662,624 (-23,196)	43,820,319 (-13,400)	26.4%



Sequence Processing

Digital Normalization

Brown et al (2014)

Retains nearly all real k -mers while discarding the majority of erroneous k -mers.

Reduces the number of over-represented k -mers (duplicates, repeats)



Sequence Processing


Uneven Sequencing Depth:

Heterozygosity


Repeated Nuclear Sequences

Organelles

Biological Contamination



Whole Genome Illumina Reads



Assembled Contigs











Zheng et al. 2013. Wheat stripe rust fungus. Nature Communications 4:2673.

Sequence Processing

Read Overlap

FLASH Majoc & Salzberg 2011 Bioinformatics 27:2957-2963.


Possible Outcomes of Paired Read Merge Algorithm

	Ground truth	FLASH Output	Label
(a)			Correct merge
(b)			Incorrect non-merge
(c)			Incorrect merge
(d)			Correct non-merge
(e)			Incorrect merge

Sequence Processing

Read Overlap


FLASH results with 1 million simulated reads, 1% error



Sequence Processing

Read Overlap

FLASH results with 1 million simulated reads, 1% error



Diverse and widespread contamination evident in the unmapped depths of high throughput sequencing data

Richard W Lusk

bioRxiv posted online January 30, 2014
Access the most recent version at doi:10.




Figure 4. Heterogeneous species appear to contaminant samples from the same tissue and experiment. The “Tumor” [18] experiment dissociated 100 individual cells from a sample of a single tumor and sequenced libraries from each. Following the analysis pipeline of a study that claimed to find different plant species in different blood plasma samples from a single experiment, I used bowtie to screen each read in each library against

Bioinformatics

Illumina data from 1 lane:

200 million paired end reads

40 gigabase

800 million lines

60 gigabyte file size

Data Transfer

FTP

HTTP

portable hard drives

[illegible]



Main Menu

HOME
CGRB BIOINFO
WEB ACCESS (ETA)
DOCUMENTATION
DONATIONS
POLICIES
ACCOUNT REQUEST
ANNOUNCEMENTS
SERVICES
WORKSHOP ACCOUNT
CLOUD STATUS
CONTACT

CGRB Bioinformatics

Resources

- > 1800 AMD and Intel processors
- >900TB shared disk space
- Hardware:
 - AMD
 - Intel
 - Dell
 - Advanced HPC
 - Hewlett Packard
 - Sun Microsystem
 - APC and Eaton
- 10G main network connection
- Gigabit private network
- Secure, climate controlled

Projects

- Assembling the Fungal Tree of Life
- Marine Microbial Genomics
- CGRB Bioinformatics
- Plantontology.org
- Phytophthora-id.org



Bioinformatics

Computing Hardware

Desktop Computer with 16 GB RAM, 1 TB HD	\$1,500
48 CPU Server with 96 GB RAM	\$9,000
24 TB Storage Array	\$10,000

Bioinformatics

Data Processing Options

Commercial Packages (Geneious, CLC Bio)

\$400 and up

Bioinformaticist or Undergraduate Programmer

NGS experience or not

Web-Services

iPlant, Galaxy



Do It Yourself

Programming (R, Python, Perl, Ruby, Java ...)

Linux + Google

Bioinformatics

Open Source Software

Blat (Kent et al 2002)
Trimmomatic (Bolger et al 2014)
khmer (Brown et al 2014)

Sequence Similarity Searching
Short Read Quality Trimming
Digital Normalization

Velvet (Zerbino & Birney 2008)
YASRA (Ratan 2009)
BWA (Li & Durbin 2009)
Abyss (Simpson et al 2009)
Trinity (Grabherr et al 2011)

Short Read Assemblers

(RNA)

SAMTools (Li et al 2009)
Tablet (Milne et al 2010)

Analyzing Mapped Reads
Visualizing Mapped Reads

Scripting Languages:

Python, R swirlstats.com

Linux

My 6 favorite commands:

sed, grep, sort, uniq, join, awk



For questions that we didn't answer:



SEQanswers

the next generation sequencing community

[SEQanswers](#) > [Bioinformatics](#)
[Bioinformatics](#)

Welcome, liston.
 You last visited: 06-27-2014 at 02:57 PM

[User CP](#)
[FAQ](#)
[Community ▼](#)
[Calendar](#)
[New Posts: 1,470](#)
[New Posts In ▼](#)
[Search ▼](#)
[Quick Links ▼](#)
[Log Out](#)

[New Thread](#)

Page 1 of 521
 [1](#)
[2](#)
[3](#)
[11](#)
[51](#)
[101](#)
[501](#)
[>](#)
[Last »](#)

Threads in Forum : Bioinformatics

[Forum Tools ▼](#)
[Search this Forum ▼](#)

	Thread / Thread Starter	Last Post ▼	Replies	Views
	Sticky: ICGC-TCGA DREAM Somatic Mutation Calling Challenge khoulahan	05-27-2014 10:47 PM by karencmartin	7	28,838
	Sticky: New Resources for 1000 Genomes (1 2) laura	12-05-2013 02:37 AM by laura	35	19,914
	Sticky: SEQwiki dan	08-09-2010 11:56 PM by scofield_gao	3	20,421
	Sticky: Software packages for next gen sequence analysis (1 2 3 ... Last Page) sci_guy	12-25-2009 06:45 PM by ECO	236	357,003
	<input checked="" type="checkbox"/> How to blast a transposon sequence to a MiSeq transposon library helelein	Today 05:06 PM by Brian Bushnell	3	84
	<input checked="" type="checkbox"/> KEGG Pathway menegidio	Today 04:21 PM by menegidio	0	27
	<input checked="" type="checkbox"/> Help with Bowtie2 aakriti	Today 04:14 PM by GenoMax	9	122
	<input checked="" type="checkbox"/> Remote blast+ skipping query sequences arundurvasula	Today 04:07 PM by GenoMax	1	69
	<input checked="" type="checkbox"/> unmark duplicates SWP	Today 03:04 PM by mebbert	2	126
	Mitochondria genome-denovo assembly bioman1	Today 02:24 PM by francicco	12	895
	<input checked="" type="checkbox"/> RNA-Seq Pathway and Gene-set Analysis Workflows in R/Bioconductor with GAGE/Pathview (1 2 3 ... Last Page) bigmw	Today 01:49 PM by tigerxu	82	7,528
	<input checked="" type="checkbox"/> BioGPS, command line? sindrle	Today 01:44 PM by sindrle	4	238
	<input checked="" type="checkbox"/> bootstrapping gene-content trees with discrete binary data using fseqboot (EBASSY) someperson	Today 01:20 PM by someperson	0	63
	<input checked="" type="checkbox"/> HMMSearch on 6-frame translated genome vs protein sequences yields different results Loddi	Today 12:51 PM by Loddi	0	54